



Toilet Safety Frame



This information is a general guideline. Please consult a healthcare provider for instructions specific to your situation.

The toilet safety frame is designed to provide support, increase stability and assist in safe toilet transfer. It is NOT designed to support the total weight of the user. Most standard toilet safety frames have a weight limit of 250#; always observe this limit.

Installation Instructions

1. Remove the toilet seat, lid, and hardware from the toilet bowl.
2. Align the holes on the mounting bracket with holes on the toilet bowl
3. Put the toilet seat and lid in place on top of the mounting bracket
4. Attach the toilet seat, lid and mounting bracket to the toilet and tighten securely.
5. On one arm support, press the snap button and insert the arm support into the mounting bracket. Slide the snap button to the inner adjustment hole. NOTE” There will be an audible “click.” Make sure the snap button fully protrudes through the hole in the mounting bracket
6. Install the other arm support in the same manner
7. To adjust the leg extension, press the snap button and slide the extension up or down to the desired hole. You will hear a “click.” Make sure both leg extensions are adjusted so the frame sits level.
8. Ensure the toilet safety frame is secure before use. Check that all snap buttons are fully protruding through the adjustment holes.